


QUICKtab Touchless Dispenser



Key Features:

- Compact, Durable Design
- Tamper-Proof TubeLOK and DispenseLOK
- Increase Productivity while Reducing Downtime and Trips for Refills
- Monitor & Track Tablet Usage for Efficient Inventory Tracking
- Lifetime Hardware Warranty



QUICKtab Sprayer



Make up a 4,306 ppm Klorese disinfectant solution by:

- ✓ Filling QUICKfoam™ with water up to the 32 oz mark.
- ✓ Add 1 - 13 gram Klorese tablet.
- ✓ Once tablet dissolves, screw on pumphead and pump blue handle to pressurize.
- ✓ Depress trigger/button to discharge Klorese foam into drain until the QUICKfoam™ level reaches the 16 oz. mark.
- ✓ Re-pressurize QUICKfoam™ as required.



QUICKfoam Drains & P-Traps



Make up a 4,306 ppm Klorese disinfectant solution by:

- ✓ Filling QUICKfoam™ with water up to the 32 oz mark.
- ✓ Add 1 - 13 gram Klorese tablet.
- ✓ Once tablet dissolves, screw on pump-head and pump blue handle to pressurize.
- ✓ Set **QUICKfoam** on counter
- ✓ Place black drain-cap over drain and press firmly to seal
- ✓ Depress trigger/button to discharge Klorese foam into drain until the QUICKfoam™ level reaches the 16 oz. mark.
- ✓ Re-pressurize QUICKfoam™ as required & move on to the next sink.



QUICKfoam-10 for Floor Drains



The Good Stuff:


- ✓ **Battery Operated**
 - ✓ Battery operated that provides 45-minutes of continuous use from one charge
- ✓ **Trigger Handle Gun**
 - ✓ Includes polypropylene trigger handle gun with a 50-degree fan tip
- ✓ **Configure With Color**
 - ✓ Hoses & lids come in a variety of colors to match your chemical management systems and quickly tells units apart



A unique concept in disinfection

Sporicidal – 4 MINUTES for *C. diff*

EPA registered Biofilm claim



- Effective Cleaner & Disinfectant across your entire facility
- Effective against ALL antibiotic resistant ESKAPE bacteria
- Safe and effective on floors
- *Cleans and disinfects in one application**
- Touchless, simple, safe and effortless
- Tested effective against Human Norovirus
- No dilution equipment to maintain or monitor
- *Sustainable*: reduced packaging & all packaging is recyclable
- Significant reduction in storage space and shipping cost
- Refill 32 oz. bottle from any water source (sink, drinking fountain, etc.**)
- "Effectively Human" – same active ingredient (HOCL) as our immune system

EPA Reg. No. 71847-7-94101



KloreSe is listed on: **EPA List N**: Disinfectants for Coronavirus (COVID-19)

EPA List K: Registered Antimicrobial for Clostridium difficile Spores

EPA List P: Registered Claims for Candida Auris

EPA Lists B, C, D, E, F, G, H, K, L, M, N, O, P, Q and Biofilm and is registered for use in electrostatic sprayers.

Remove & Eliminate: PerfectCLEAN ✓ & KLORESE - the perfect partnership

* *C. diff, Ca. auris and OSHA BBP requires pre clean*

** *Effective in water up to 400 pm hardness*

Item No.	Description	EPA Reg. No. 71847-7-94101	Packaging
Klorese-110	110 - Foil Wrapped 13.1 g Tablets – makes 110 quarts		2 PAILS per CASE
Klorese-960	48 - QUICKtab Dispenser Tubes in 4 cases – makes 960 Quarts		48 - 20 TABLET TUBES per CASE
K632	6 - Pack 32 oz pre-measured Bottles with Flip Top for 13.1 g Tablets		6 BOTTLES per CASE
K6SPRAY	6 - Pack Green/White		6 SPRAYERS per CASE
K25LAB	25 - Adhesive Backed "6 Day Expiration" Labels for 32 oz Re-Fillable Bottle		25 LABELS per POLYBAG